

AED Plus - Volautomatisch - Franstalig (incl. batterij en elektroden)

08/10/2025 00:54:57



Artikelnummer	3100FR
Prijs excl. BTW	€ 1.575,69
Prijs incl. BTW	€ 1.906,58
Voorraad	Nog slechts 5 op voorraad

Afbeeldingen



	 ZOLL AED Plus [®]	 POWERLIFE 50	 ZOLL AED 3P
Model overview			
Detection of a RCP or arrest and confirmation of successful resuscitation	Real CPR help [®]	Estimation of a RCP coefficient 0-100	Real CPR help [®] available
Response to profound low impedance	+	+	+
Instructions: sounds and text-to-speech	+	+	+
Images of scenes / graphics	Screens: graphics		Images on screen
Screen: touch			+
Monitor of rhythm via RCP			+
Response time	12 seconds	12 seconds	8 seconds
Maximum pathology	+ Algorithms of analysis and diagnosis available for 17 values + Detection of pathological values in	+ Recognition of rhythm in case of ventricular fibrillation, pulseless tachycardia + Four other pathological rhythms	+ Algorithms of analysis and diagnosis available for 17 values + Universal CPR feedback "per minute and value"
Other features	+		
Key features			
Duration of use in life	1200 (typical) - 3 ans	Overcharge intelligent [®] - 3 ans	Overcharge CRM intelligent [®] - 3 ans
Depth of use in life of the battery	3 ans	Battery intelligent de 3 ans	Battery intelligent de 3 ans
Indicator of charge level	+	+	+
Accessories and others	Indication module	Overcharge	Indication module
Compatible with the Wi-Fi			+
Connectivity and data management and storage	WiFi	WiFi	Wi-Fi or USB
Security	+ IP50	+ IP54 + Test de chute à 1 mètre	+ IP54 + Test de chute à 1 mètre
Options for programming	+	+	+